

SRM
INSTITUTE OF SCIENCE & TECHNOLOGY
(Chartered as the University with effect from 1st Aug, 2010)

CHEMISTRY


44th and 45th class, 11-11-2021

Dr. K. Ananthanarayanan
Associate Professor (Research)
Department of Chemistry
Room No 319, 3rd Floor, Research Building

Email : ananthak@srmist.edu.in
Phone : 9840154665

11/11/2021

Page 1



SRM
INSTITUTE OF SCIENCE & TECHNOLOGY
(Chartered as the University with effect from 1st Aug, 2010)

Virtual lab.

Determination of strength of an acid
using pH meter

Expt. No. : 7

11/11/2021

2

Experiment



□ Aim :

To find out the strength of given hydrochloric acid solution by titrating it against sodium hydroxide (0.05N) using pH meter

□ Materials required:

pH meter, electrodes, beaker, 100mL standard flask, pipette, burette, funnel, glass rod

□ Chemicals required :

Hydrochloric acid, NaOH solution, distilled water

11/11/2021

3

pH of a solution




□ The pH of a solution is the negative logarithm of the hydrogen ion concentration.

□ pH may be represented mathematically as :

$$\text{pH} = -\log_{10} [\text{H}^+]$$

□ If the $[\text{H}^+]$ concentration is greater than $1 \times 10^{-7} \text{ M}$, pH is less than 7 and vice-versa.

□ $\text{pH} < 7$ indicate **acidic solution**, $\text{pH} > 7$ indicate **basic solution** and $\text{pH} = 7$ means **neutral solution**.



11/11/2021

4

Potentiometric titration



- ☐ **Volumetric method** for the measurement of **changes in the e.m.f between two electrodes on the addition of titrant of known concentration to a solution of unknown concentration.**
- ☐ In this potentiometric titration setup an **indicator electrode (ISE electrode) and reference electrode (calomel) is coupled to form electrochemical cell.**
- ☐ The overall cell potential is calculated as :

$$E_{\text{cell}} = E_{\text{ind}} - E_{\text{ref}} + E_{\text{sol}}$$

- ☐ E_{sol} is the **potential drop over the test solution** between the two electrodes which is recorded at intervals as the titrant is added.

11/11/2021

5

Potentiometric titration



- ☐ We have seen that the potential of a metallic electrode arises from the tendency of an oxidation-reduction reaction to occur at an electrode surface (redox titration).
- ☐ In membrane electrodes, in contrast, the observed potential is a kind of junction potential that develops across a membrane that separates the analyte solution from a reference solution.
- ☐ **Requirements:** minimal solubility, need some electrical conductivity and selective reactivity with the analyte

11/11/2021

6

pH meter



11/11/2021


7

pH meter



❑ A pH meter consists of two electrodes: **a glass electrode (ISE)**, which is sensitive to the concentration of hydrogen ions in solution, and **a calomel reference electrode**.

❑ The reference electrode supplies a constant potential ($E^\circ = +0.24 \text{ V}$ versus the SHE) as determined by the half-reaction :




❑ The pH meter measures the total potential across the two electrodes and displays this measurement on a scale calibrated in pH units

11/11/2021

8

Principle



11/11/2021

9

Principle



- ❑ In order to measure the voltage, two points with different electrical potential values are required.
- ❑ The reference electrode is designed to maintain a constant electrical potential that is independent of the sample composition and temperature.
- ❑ In contrast, the hydrogen ion selective electrode (ISE) with glass membrane provides an electrical potential that is dependent upon the activity of hydrogen (H^+) ions in the sample solution.
- ❑ Therefore, both the reference electrode and the hydrogen ISE are needed when determining pH.


11/11/2021

10

Principle




- ❑ This membrane is filled with a buffer solution of known pH (typically pH = 7).
- ❑ This electrode design creates an environment with constant binding of H^+ ions on the inside of the glass membrane, while the outside of the glass membrane is exposed to the sample where a variable amount of H^+ ions exist.
- ❑ The difference in H^+ ions creates a potential that is read versus the stable potential of the reference electrode.



11/11/2021

11


Principle



11/11/2021

12


Principle



● Si ○ Oxygen ○ Cations

Structure of Silica pH Sensitive Glass


The ideal mixture was found to be 72% silicon dioxide, 6% calcium oxide and 22% disodium oxide.



11/11/2021

13

Principle



11/11/2021

14

Principle



- ❑ The glass used in a pH electrode contains sodium (Na^+) ions.
- ❑ When placed in contact with a solution, a microscopically thin (10 nm thick) hydrated gel layer forms on the surface of the glass.
- ❑ Within this gel layer, hydrogen ions from the solution can migrate a very short distance into the glass and displace sodium ions.
- ❑ The number of H^+ ions that migrate into the gel layer depends directly on the concentration of H^+ in the solution.

11/11/2021

15

Principle



- ❑ For the **internal solution, the concentration of H^+ is constant** so a fixed number of Na^+ ions are displaced from the gel layer by H^+ ions.
- ❑ For the **external solution, the number of H^+ ions changes with the pH** so the number of Na^+ ions displaced from the gel layer varies as a function of the pH.
- ❑ The ratio of **$[\text{H}^+]$ to $[\text{Na}^+]$ in the outside gel layer influences the magnitude of the membrane potential.**
- ❑ The reason the **membrane potential varies** with the ratio of the two ions is because the **smaller H^+ ions have a higher mobility than the larger Na^+ ions.**

11/11/2021

16

Principle



- A voltmeter measures the voltage generated by the solution and displays it as a pH-measurement.
- An increase in voltage means more hydrogen ions and an increase in acidity, so the meter shows it as a decrease in pH; in the same way, a decrease in voltage means fewer hydrogen ions, more hydroxide ions, a decrease in acidity, an increase in alkalinity, and an increase in pH.
- \uparrow voltage = more H^+ /less OH^- = \uparrow acidity = \downarrow pH
- \downarrow voltage = less H^+ /more OH^- = \downarrow acidity = \uparrow pH

11/11/2021

17

Principle



- In order to measure pH, the cell is set up –

Ag | AgCl, 0.1.M HCl | glass | Unknown solution | sat.KCl | Hg₂Cl₂ | Hg

The glass electrode and calomel electrode are dipped in a solution of known pH.

$$E_{\text{cell}} = E_{\text{ref(int)}} - E_{\text{ref(samp)}} + E_{\text{mem}} + E_j$$

11/11/2021

18

Principle



$$E_{\text{mem}} = E_{\text{asym}} - \frac{RT}{zF} \ln \frac{(a_A)_{\text{int}}}{(a_A)_{\text{samp}}}$$

$$0.0592 \log [(A_{\text{inner}}) / (A_{\text{outer}})]$$

$(a_A)_{\text{samp}}$ is the analyte's concentration in the sample, $(a_A)_{\text{int}}$ is the concentration of analyte in the ion-selective electrode's internal solution, and z is the analyte's charge. Ideally, E_{mem} is zero when $(a_A)_{\text{int}} = (a_A)_{\text{samp}}$. The term E_{asym} , which is an asymmetry potential, accounts for the fact that E_{mem} is usually not zero under these conditions.

11/11/2021

19

Principle



- ☐ The **pH values of the solution** at different stage of acid–base neutralization are determined and **plotted against the volume of alkali added.**
- ☐ On **adding a base to an acid**, the **pH rises slowly** in the initial stages as the concentration of H^+ ion decreases gradually.
- ☐ At the **equivalence point, it increases rapidly** as at the equivalent point H^+ ion concentration is very small and then it flattens out after the end point.
- ☐ The end point of the titration can be detected where the pH value changes most rapidly.

11/11/2021

20

Procedure



- ☐ Standardize the pH meter using different buffers of known pH, wash the electrode with distilled water.
- ☐ The given acid is diluted to 100 ml using distilled water. **10 ml** of this made up solution is pipetted out into a 250 mL clean beaker.
- ☐ 100 ml of distilled water is added to it, so that the electrode is completely dipped. Note the initial pH of the pure acid solution.
- ☐ Fill the burette with standard NaOH solution (**0.05 N**) and run down into the beaker in small increments (1 mL). Stir the solution well using glass rod.

11/11/2021

21

Procedure contd..



- ☐ Note down the pH after every successive addition.
- ☐ Continue the titration till beyond the neutralization point as indicated by an abrupt change in pH. (at least 15-20 increments).
- ☐ Plot a pilot graph of volume of NaOH versus pH. The midpoint of the **'S' shaped curve** of the graph gives the equivalence point.
- ☐ In order to get an accurate end point, perform one more similar titration in a small volume range (3.0 mL on either side of the abrupt change in pH)


11/11/2021

22


Procedure contd..



- ☐ Measure the pH after addition of every 0.2 mL of standard NaOH solution.
- ☐ Plot a fair graph of volume of NaOH versus $\Delta \text{pH} / \Delta V$.
- ☐ Find out the exact end point from the fair graph. The peak point of the curve from the fair graph gives the end point



Pilot Graph



Fair Graph

11/11/2021

23

Tabular column



Table-1 : Pilot titration (HCl Vs NaOH)

S.No	Volume of NaOH (V) ml	pH
1.		
2.		
3.		
....		
16		

11/11/2021

24

Tabular column



Table-2 : Fair titration (HCl Vs NaOH)

S.No.	Volume of NaOH (V) ml	pH	Δ pH	ΔV	$\Delta \text{pH} / \Delta V$
1.					
2.					
3.					
....					
20.					

11/11/2021

25

Calculation



Volume of NaOH V_1 = ----- ml [from FAIR graph]

Strength of NaOH N_1 = ----- N

Volume of HCl V_2 = 10 ml

Strength of HCl N_2 = ----- N

Therefore strength of HCl (N_2) = $V_1 \times N_1 / V_2$

= -----N.

11/11/2021

26

Result



□ The strength of the given HCl solution = N

11/11/2021

27

Table 1 Determination of the strength of an acid using pH meter

Sl.No.	Volume of NaOH (mL)	pH
1	0	2.84
2	1	2.86
3	2	2.90
4	3	2.94
5	4	3.00
6	5	3.07
7	6	3.15
8	7	3.26
9	8	3.40
10	9	3.61
11	10	4.13
12	11	7.11
13	12	9.85
14	13	10.41
15	14	10.67
16	15	10.84



11/11/2021

28

Table 2 Determination of the strength of an acid using pH meter

Sl.No	Volume of NaOH (mL)	pH	Δ pH	Δ V (mL)	Δ pH/ Δ V
1	9	3.69	-----	-----	
2	9.2	3.73	0.04	0.2	0.2
3	9.4	3.83			
4	9.6	3.91			
5	9.8	4.07			
6	10	4.31			
7	10.2	4.86			
8	10.4	5.83			
9	10.6	6.29			
10	10.8	6.80			
11	11	7.20			
12	11.2	8.70			
13	11.4	9.26			
14	11.6	9.47			
15	11.8	9.72			
16	12	9.87			



29



17	12.2	10.05			
18	12.4	10.17			
19	12.6	10.25			
20	12.8	10.33			
21	13	10.42			
22	13.2	10.54			
23	13.4	10.72			

11/11/2021

30