



CHEMISTRY

50th and 51st class Date : 18-11-2021

Dr. K. Ananthanarayanan
Associate Professor (Research)
Department of Chemistry
Room No 319, 3rd Floor, Raman Research Park

Email : ananthak@srmist.edu.in

Phone : +91-9840154665



In this class...

Determination of molecular weight of a polymer
by viscosity average method

Expt. No. : 8

Experiment



□ Aim :

To determine the molecular weight of a polymer in solution by using a Ostwald viscometer.

□ Materials required:

Ostwald viscometer, beaker, stop watch, standard flasks, pipette, suction bulb

□ Chemicals required :

Given polymer solution, suitable solvents


11/18/2021

3

Polymer



- A polymer is a substance or material consisting of very large molecules, or macromolecules, composed of many repeating subunits.
- Polymerization is the process by which monomers (smaller chemical units) are combined to form a polymer.





11/18/2021

2nd paper presentation-Vidya

4

Polymer used – polyvinyl alcohol



11/18/2021

5

Molecular weight



- ☐ The molecular weight of a synthetic polymer **does not have a single value, since different chains will have different lengths and different numbers of side branches.**
- ☐ There will **therefore be a distribution of molecular weights, so it is common to calculate the average molecular weight** of the polymer.
- ☐ Several different ways to define the average molecular weight - the **two most common being the number average molecular weight and the weight average molecular weight.**
- ☐ Other averages exist, such as the **viscosity average molecular weight.**

11/18/2021

6

Number average molecular weight

The number average molecular weight **is defined as the total weight of polymer divided by the total number of molecules.**

$$\text{Total weight} = \sum_{i=1}^{\infty} N_i M_i$$

where N_i is the number of molecules with weight M_i

$$\text{Total number} = \sum_{i=1}^{\infty} N_i$$

The number average molecular weight is therefore given by:

$$\overline{M}_N = \frac{\sum_{i=1}^{\infty} N_i M_i}{\sum_{i=1}^{\infty} N_i}$$

11/18/2021

7

Weight average molecular weight

The weight average molecular weight **depends not only on the number of molecules present, but also on the weight of each molecule.** To calculate this, N_i is replaced with $N_i M_i$.

$$\overline{M}_W = \frac{\sum_{i=1}^{\infty} N_i M_i^2}{\sum_{i=1}^{\infty} N_i M_i}$$

This can also be written as:

$$\overline{M}_W = \sum_{i=1}^{\infty} w_i M_i$$

where w_i is the weight fraction of polymer with molecular weight M_i .

11/18/2021

8

Molecular weight



Consider a polymer sample comprising of 5 moles of polymer molecules having molecular weight of 50,000 g/mol and 9 moles of polymer molecules having molecular weight of 30,000 g/mol.

$$\bar{M}_n = \frac{(9 \text{ mol} \times 30,000 \text{ g/mol}) + (5 \text{ mol} \times 50,000 \text{ g/mol})}{9 \text{ mol} + 5 \text{ mol}} = 37,000 \text{ g/mol}$$

$$\bar{M}_w = \frac{9 \text{ mol}(30,000 \text{ g/mol})^2 + 5 \text{ mol}(50,000 \text{ g/mol})^2}{9 \text{ mol}(30,000 \text{ g/mol}) + 5 \text{ mol}(50,000 \text{ g/mol})} = 40,000 \text{ g/mol}$$

Viscosity average molecular weight



$$\bar{M}_V = \left[\frac{\sum_i N_i M_i^{1+a}}{\sum_i N_i M_i} \right]^{1/a}$$

where a is the exponent in Staudinger-Mark-Houwinks equation: $[\eta] = k \cdot M^a$ ($[\eta]$ = intrinsic viscosity). For a flexible polymer in a good solvent, $a \approx 0.7$.

\bar{M}_V is determined by viscosity measurements.

Determination of molecular weight



a. Absolute method :

Mass spectrometry
Colligative property
End group analysis
Light scattering
Ultracentrifugation.


b. Relative method : Solution viscosity

c. Fractionation method : GPC (Gel Permeation Chromatography)


Viscosity



- ☐ Viscosity is a measure of a fluid's resistance to flow.
- ☐ It describes the internal friction of a moving fluid caused due to intermolecular forces between particles within the fluid.
- ☐ A fluid with large viscosity resists motion because its molecular makeup gives it a lot of internal friction.



Laminar flow



- ❑ Laminar flow occurs in layers without mixing. Notice that viscosity causes drag between layers as well as with the fixed surface.
- ❑ An obstruction in the vessel produces turbulence. Turbulent flow mixes the fluid. There is more interaction, greater heating, and more resistance than in laminar flow.


11/18/2021

13

Poiseuille's equation




- ❑ The precise definition of viscosity is based on **laminar**, or non turbulent, flow characterized by the smooth flow of the fluid in layers that do not mix.



- ❑ The **flow rate Q** of the laminar flow of an incompressible fluid having **viscosity \$\eta\$** through a horizontal tube of **uniform radius r** and **length l** is given by the Poiseuille's equation :

$$Q = \frac{(P_2 - P_1)\pi r^4}{8\eta l}$$



11/18/2021

14

Principle



- In a capillary viscometer, the viscosity of a liquid is proportional to the time taken by a known volume of liquid to flow through the capillary under a specified hydrostatic pressure at a fixed temperature.
- Provided the flow is laminar, using the Poiseuille's equation it is possible to show that if t , η and ρ are the flow time, viscosity and density of solution :

$$\frac{\eta}{\eta_o} = \frac{\rho}{\rho_o} \frac{t}{t_o}$$

$$\frac{\eta}{\eta_o} = \text{Relative viscosity } \eta_{rel}$$

The specific viscosity is defined as $\eta_{sp} = \eta_{rel} - 1$

11/18/2021

15



Relative viscosity : $\eta_{rel} = \frac{\eta}{\eta_o} = \frac{t}{t_o}$

η : solution viscosity
 η_o : solvent viscosity
 t : flow time of solution
 t_o : flow time of solvent


Specific viscosity : $\eta_{sp} = \frac{\eta - \eta_o}{\eta_o} = \frac{t - t_o}{t_o} = \eta_{rel} - 1$

Reduced viscosity : $\eta_{rel} = \frac{\eta_{sp}}{c}$

Intrinsic viscosity : $[\eta] = \lim_{c \rightarrow 0} \frac{\eta_{sp}}{c}$

- The ratio of a solution's specific viscosity to the concentration of the solute, extrapolated to zero concentration.
- Intrinsic viscosity reflects the capability of a polymer in solution to enhance the viscosity of the solution

Ostwald Viscometer




Principle contd..

- ❑ The Staudinger – Mark – Houwink's equation which relates η_{int} with molecular weight

$$\eta_{int} = K(M)^\alpha$$

- ❑ Where 'K' is an empirical parameter characteristic of a particular solute – solvent pair and 'α' is a shape parameter. (both are constants). a = Mark-Houwink constant relating to structure (0 to 0.1 – sphere, 0.35 to 0.80 – random coil, 1.5 to 2 – rigid rod).
- ❑ Plot of η_{sp}/C Vs C , extrapolated to C= 0 gives η_{int}
- ❑ From known values of K and α , molecular weight can be determined.



Procedure



- ☐ 1% solution of polymer (PVA) in water will be supplied. Prepare at least 5 dilutions i.e 0.1%, 0.2% , 0.3% , 0.4% and 0.5% polymer in water.
- ☐ Dilutions can be done by using the volumetric expression $N_1V_1 = N_2V_2$.
- ☐ Set up the Ostwald viscometer and measure the flow time (t_0) of a fixed volume of pure solvent (12.5 ml). Take the average reading (at least three iterations).
- ☐ Rinse the viscometer thoroughly with the most dilute solution, measure the flow time (t_1) keeping the flow volume the same.
- ☐ Repeat the procedure for other solutions.
- ☐ Calculate η_{rel} and η_{sp} . Plot η_{sp}/C Vs C , extrapolate to $C=0$ to obtain η_{int} . From the given values of K and α , calculate the molecular weight.

11/18/2021

19

Tabular column



S.No.	Concentration of the polymer solution	Time of flow in sec (average)	Relative viscosity $\eta/\eta_0 = t/t_0$	Specific viscosity $\eta_{sp} = \eta/\eta_0 - 1$	Reduced viscosity η_{sp}/C
1.	Pure Solvent	t_0			
2.	0.1%	t_s			
3.	0.2%	t_s			
4.	0.3%	t_s			
5.	0.4%	t_s			
6.	0.5%	t_s			

11/18/2021

20

Calculation



Solvent used Water


$$\eta_{\text{int}} = K \times M^{\alpha}$$

$$\log \eta_{\text{int}} = \log K + \alpha \log M$$

$$\alpha \log M = \log \eta_{\text{int}} - \log K$$

$$\log M = \frac{\log \eta_{\text{int}} - \log K}{\alpha}$$

$$M = \text{Anti log} \frac{[\log \eta_{\text{int}} - \log K]}{\alpha}$$



Plot η_{sp}/C Vs C .

Extrapolate the straight line to Y axis

The intercept is η_{int} .

From the relationship, $\eta_{\text{int}} = KM^{\alpha}$ where K and α are constants,

(For PVA solution $K = 45.3 \times 10^{-3}$, $\alpha = 0.64$)

11/18/2021

21

Result



□ The molecular weight of the given polymer is =

11/18/2021

22

Tabular column, values



Sl. No	Concentration of the polymer solution (c) (%)	Time of flow in sec (average)	Relative viscosity $\frac{\eta}{\eta_0} = \frac{t}{t_0}$	Specific viscosity $\eta_{sp} = (\text{relative viscosity} - 1)$	Reduced viscosity $\frac{\eta_{sp}}{c}$
1	Pure solvent	38 (t_0)	-----	-----	-----
2	0.1	43	1.131	0.131	1.31
3	0.2	46			
4	0.3	50			
5	0.4	55			
6	0.5	62			