



CHEMISTRY

16th and 17th class, Date : 04-10-2021

Dr. K. Ananthanarayanan
Associate Professor (Research)
Department of Chemistry
Room No 319, 3rd Floor, Raman Research Park

Email : ananthak@srmist.edu.in

Phone : 9840154665



In this class..

Estimation of total hardness, permanent
and temporary hardness by EDTA
method

Expt. No. : 2

Experiment



□ Aim :

To estimate the amount of total hardness, permanent hardness and temporary hardness of a given sample of water by EDTA method using ammonia buffer (pH =10) and Eriochrome black-T indicator.

□ Materials required:

Pipette, burette, conical flask, standard flask, funnel, beaker.

Chemicals required :

EDTA solution, standard hard water, sample water, Eriochrome black T indicator (EBT), $\text{NH}_3\text{-NH}_4\text{Cl}$ buffer solution (pH 10).

5 October 2021

Dr K Ananthanarayanan

3

Hardness of water



□ Hard water is water that contains cations with a charge of +2, especially Ca^{2+} and Mg^{2+}

□ At a **level over 100 mg/L** measured as calcium carbonate, a given water source is considered hard. (Not potable)

Hardness (mg/L)	Degree of hardness
0-75	Soft
75-100	Moderately hard
150-300	Hard
>300	Very hard

5 October 2021

Dr K Ananthanarayanan

4

Hardness of water, scale formation



5 October 2021

Dr K Ananthanarayanan

5

Types of hardness



☐ **Temporary hardness:**

It is due to the presence of bicarbonates of calcium and magnesium and can be removed by boiling.

☐ **Permanent hardness:**

It is due to the presence of sulphates, nitrates and chlorides of calcium and magnesium.

☐ Hardness is usually reported as **parts per million (ppm) of calcium carbonate (by weight).**

☐ A water supply with a hardness of **100 ppm contains** the equivalent of **100 g of CaCO_3 in 1 million g of water or 0.1 g in 1 L of water**

5 October 2021

Dr K Ananthanarayanan

6

Hardness, representation



Clark's degree ($^{\circ}\text{Cl}$): The number of parts of calcium carbonate equivalent hardness present in **70,000 or (7×10^4)** parts of water.

1 $^{\circ}$ Clark = 1 part of CaCO_3 eq. hardness per 70,000 parts of water. (**14.254 ppm**)

Degree French ($^{\circ}\text{Fr}$): The number of parts of calcium carbonate equivalent hardness presents in **10^5 parts of water**.

1 $^{\circ}$ Fr = 1 part of CaCO_3 hardness eq per 10^5 parts of water. (**10 ppm**)

5 October 2021


Dr K Ananthanarayanan

7

Complexometric titration



- ☐ A titration based on the formation of coordination complexes between a **metal ion and complexing agent** (or chelating agent) to form soluble complexes. (**Hardness in water**)
- ☐ Complex-forming reactions involving **many metal ions** can serve as a basis for accurate and convenient titrations for such metal ions. **High accuracies** and offer the possibility of determinations of metal ions at the **millimole levels**.



5 October 2021



Dr K Ananthanarayanan

8

Principle



- The hardness causing metal ions form a **wine-red coloured weak complex with Eriochrome black-T indicator** in the presence of buffer solution.
- When EDTA is added, the **indicator is replaced by EDTA and a stable complex is formed**. Due to **liberation of Eriochrome black-T indicator**, wine **red colour changes to steel blue** which marks the end point of titration.



5 October 2021


Dr K Ananthanarayanan

9

Role of Eriochrome black-T



- Eriochrome black T is an azo dye used as a complexometric indicator.
- In its **deprotonated form, Eriochrome black-T is blue**. It turns **red when it forms a complex** with calcium, magnesium, or other metal ions.



- The **characteristic blue end-point is reached** when sufficient EDTA is added and the metal ions bound to the indicator are chelated by EDTA, leaving the free indicator molecule.

5 October 2021




Dr K Ananthanarayanan

10

Role of EDTA



- EDTA is a **hexadentate complexing or chelating agent** used to capture the metal ions.




- Capturing metal ions causes water to become softened, but metal ions are not completely removed from water.
- EDTA simply binds the metal ions to it very tightly thus forming a **strong and stable metal complex**.

5 October 2021

Dr K Ananthanarayanan

11

Role of EDTA




EDTA (anionic form)

5 October 2021

Dr K Ananthanarayanan

12

Complexometric titration



- ❑ The **Y⁴⁻ ion that forms a 1:1 complex** with the metal ion is the completely deprotonated anion of ethylenediaminetetraacetic acid ("**H₄Y**").
- ❑ At **pH=10**, the EDTA is present in solution primarily as its monoprotonated form, **HY³⁻**.

5 October 2021


Dr K Ananthanarayanan

13


Role of ammonia buffer



- ❑ The buffer solution is used to resist the change in pH as all reactions between **metal ions and EDTA are pH dependent**.
- ❑ With increasing the pH, each hydrogen ion in the carboxyl groups of EDTA will start to dissociate. Above pH 10, EDTA⁴⁻ is predominant. As we need EDTA⁴⁻ to react with the metal ions present in the titration solution, we use **pH 10 buffer** such as ammonia/ammonium chloride.



EDTA 4-




EDTA - metal complex

5 October 2021

Dr K Ananthanarayanan

14

Role of pH




- EDTA versus metal ion titration, if $pM (= -\log[M^{n+}])$, where M^{n+} signifies the metal ion whose concentration is required) is plotted against the volume of EDTA solution added

5 October 2021

Dr K Ananthanarayanan

15


Formation constant and equilibrium



with a corresponding equilibrium constant for **formation** expressed as

$$K_f (Ca-EDTA^{2-}) = \frac{[Ca-EDTA^{2-}]}{[Ca^{2+}][EDTA^{4-}]}$$

Formation constant is very large 2.9×10^{16} ,
its equilibrium position lies far to the right.





5 October 2021


Dr K Ananthanarayanan

16

Equation, calculation



So if we let M^{2+} be the total of Ca^{2+} and Mg^{2+} ions in solution then



Note that 1 mole of M^{2+} reacts with 1 mole of $\text{H}_2\text{EDTA}^{2-}$.
At the equivalence point of the titration:

$$\begin{aligned} \text{moles } \text{M}^{2+} &= (\text{moles } \text{Ca}^{2+} + \text{moles } \text{Mg}^{2+}) \\ &= \text{moles } \text{H}_2\text{EDTA}^{2-} \end{aligned}$$

5 October 2021

Dr K Ananthanarayanan

17

CaCO_3 as reference for hardness



- ☐ CaCO_3 is stable, non-hygroscopic and is obtained in pure form. Therefore a standard hard water solution can be prepared by dissolving accurately weighed CaCO_3 in HCl and can make up to a known volume.
- ☐ CaCO_3 is insoluble in water. Therefore it can be easily precipitated in water treatments.
- ☐ Molecular weight of CaCO_3 is 100, so mathematical calculations are easy.

1.0 gm of pure CaCO_3 dissolved in minimum quantity of conc. HCl and diluted to a one litre with distilled water.

5 October 2021

Dr K Ananthanarayanan

18

Procedure



Standardisation of EDTA

- ☐ Pipette out 20mL of standard hard water into a clean conical flask.
- ☐ Add 5mL of the buffer solution and 3 or 4 drops of eriochrome black-T indicator. The solution turns wine red in colour.
- ☐ Titrate the wine red coloured solution against EDTA taken in the burette.
- ☐ The change from wine red to steel blue colour is the end point. Repeat the titration for concordant values. Let the titer value be V_1 mL

5 October 2021

Dr K Ananthanarayanan

19

Procedure



Determination of total hardness

- ☐ Pipette out 20mL of sample hard water into a clean conical flask.
- ☐ Add 5mL of the buffer solution and 3 or 4 drops of eriochrome black-T indicator. The solution turns wine red in colour.
- ☐ Titrate the wine red coloured solution against EDTA taken in the burette.
- ☐ The change from wine red to steel blue colour is the end point.
- ☐ Repeat the titration for concordant values. Let V_2 mL be the volume of EDTA consumed.

5 October 2021

Dr K Ananthanarayanan

20

Procedure



Determination of permanent hardness

- ☐ Take 100 mL of hard water sample in a 250 mL beaker and boil gently for about one hour. (**boiled water will be given**)
- ☐ Cool and filter into a 100mL standard flask and make the volume up to the mark. Take **20mL of this solution** and proceed the titration in the same way.
- ☐ The volume of EDTA consumed corresponds to the permanent hardness of the water sample. Let the titer value be **V₃ mL.**
- ☐ Temporary hardness is calculated by subtracting **permanent hardness from total hardness.**

5 October 2021

Dr K Ananthanarayanan

21

Tabular column



Standardisation of EDTA

Table 1
Std Hard water Vs EDTA solution

S.No.	Vol. of Std. hard water ml	Burette Reading		Vol. of EDTA ml	Indicator
		Initial	Final		
					EBT

5 October 2021

Dr K Ananthanarayanan

22

Calculation, standardisation



1 ml of Standard hard water	=	1 mg of CaCO_3
Volume of standard hard water taken	=	20ml
20ml of Standard hard water	=	20mg of CaCO_3
Volume of EDTA consumed	=	V_1 ml (from table-1)
V_1 ml EDTA solution	=	20mg CaCO_3
Therefore 1ml EDTA will be	=	$\frac{20}{V_1}$ mg of equivalent CaCO_3
		A mg CaCO_3

5 October 2021

Dr K Ananthanarayanan

23

Tabular column



Determination of total hardness

Table 2
Sample Hard water Vs EDTA Solution

S.No.	Vol. of Sample hard water ml	Burette Reading		Vol. of EDTA ml	Indicator
		Initial	Final		
					EBT

5 October 2021

Dr K Ananthanarayanan

24

Determination of Total hardness

Volume of EDTA consumed = **X mL (V_2 from Table-2)**

1 mL of EDTA = $20/V_1$ [mg CaCO_3] = **A mg CaCO_3**

V_2 mL of EDTA = **A $\times V_2$ [mg CaCO_3] = B mg of CaCO_3**

If 20 mL of sample hard water taken for titration = **B mg of CaCO_3**

Then, 1000 mL will contain = **(B \times 1000 mg CaCO_3) / 20**

Total hardness = TH ppm

Tabular column

Determination of permanent hardness

Table 3
Boiled sample hard water Vs EDTA solution

S.No.	Vol. of boiled water ml	Burette Reading		Vol. of EDTA ml	Indicator
		Initial	Final		
					EBT



Determination of Permanent hardness:

Volume of EDTA consumed = **X mL (V_3 from Table-3)**

1 mL of EDTA = $20/V_1$ [mg CaCO_3] = A mg CaCO_3

V_3 mL of EDTA = $A \times V_3$ [mg CaCO_3] = **C mg of CaCO_3**

If 20 mL of sample hard water taken for titration = **C mg of CaCO_3**

Then, 1000 mL will contain = $(C \times 1000 \text{ mg } \text{CaCO}_3) / 20$

= PH ppm

Permanent hardness = ppm

5 October 2021

Dr K Ananthanarayanan

27



Estimation of Temporary hardness

Temporary Hardness = Total hardness – Permanent hardness

RESULT:

☐ The total hardness of sample hard water is =

☐ The permanent hardness of sample hard water is =

☐ The temporary hardness of sample hard water is =

5 October 2021

Dr K Ananthanarayanan

28

Limitations of EDTA titration

1. Many EDTA titrations are carried out under alkaline pH which may lead to **formation of insoluble hydroxides** or basic salts that may compete with the complexation process.
2. Since EDTA forms stable complexes with most of the metal ions, **it lacks selectivity** if it is used to estimate a single metal cations from a solution of mixture of metal ions.
3. The change in equilibrium from metal-indicator complex to the metal-EDTA complex should be **sharp and rapid**.
4. Conditional formation constant of the metal-EDTA complex K'_{MY} to K'_{In} should **be of the order 10^4** to provide a good end point

5 October 2021

Dr K Ananthanarayanan

29

Table-1 Standard Hard water vs EDTA Solution

Sl.No.	Volume of Sample Hard water (mL)	Burette Reading (mL)		Volume of EDTA solution (mL)	Indicator
		Initial	Final		
1	20	0	22.9	22.9	Eriochrome Black-T
2	20	0	22.9		

5 October 2021

Dr K Ananthanarayanan

30

Table-2 Sample Hard water vs EDTA Solution

Sl. No.	Volume of Sample Hard water (mL)	Burette Reading (mL)		Volume of EDTA solution (mL)	Indicator
		Initial	Final		
1	20	0	10.8	10.8	Eriochrome Black-T
2	20	0	10.8		

5 October 2021

Dr K Ananthanarayanan

31

Table-3 Boiled sample water vs EDTA Solution

Sl.No.	Volume of Sample Hard water (mL)	Burette Reading (mL)		Volume of EDTA solution (mL)	Indicator
		Initial	Final		
1	20	0	7.1	7.1	Eriochrome Black-T
2	20	0	7.1		

5 October 2021

Dr K Ananthanarayanan

32



Determination of Total hardness:

Volume of EDTA consumed = 10.8 mL (V_2 from Table-2)

$$1 \text{ mL of EDTA} = 20/V1 \text{ mg CaCO}_3 = \mathbf{0.873 \text{ mg CaCO}_3}$$
$$V_2 \text{ mL of EDTA} = \frac{20}{V_1} \times V_2 \text{ mg CaCO}_3 = 9.432 \text{ mg of CaCO}_3$$

If 20 mL of sample hard water taken for titration =

$$20/V_1 \times V_2 \text{ mg CaCO}_3$$

Then, 1000 mL will contain

$$= (20/V_1 \times V_2 \times 1000) / 20$$
$$= 0.471.6 \times 1000$$

Total hardness = 471.61 ppm

5 October 2021

Dr K Ananthanarayanan

33



Determination of Permanent hardness:

Volume of EDTA consumed = 7.1 mL (V_3 from Table-3)

$$1 \text{ mL of EDTA} = 20/V_1 \text{ mg CaCO}_3 = 0.873 \text{ mg CaCO}_3$$
$$V_3 \text{ mL of EDTA} = 20/V_1 \times V_3 \text{ mg CaCO}_3 = \mathbf{6.198 \text{ mg of CaCO}_3}$$

If 20 mL of sample hard water taken for titration

$$= 20/V_1 \times V_3 \text{ mg CaCO}_3$$

Then, 1000 mL will contain = $(20/V_1 \times V_3 \times 1000) / 20$

$$= 309.91 \text{ ppm}$$

Permanent hardness = 309.91 ppm

5 October 2021

Dr K Ananthanarayanan

34

Estimation of Temporary hardness

$$\begin{aligned}\text{Temporary Hardness} &= \text{Total hardness} - \text{Permanent hardness} \\ &= 471.60 \text{ ppm} - 309.91 \text{ ppm} \\ &= 161.69 \text{ ppm}\end{aligned}$$

Result:

Total Hardness = 471.61 ppm

Permanent Hardness = 309.91 ppm

Temporary Hardness = 161.70 ppm

Thank you all for your attention

Information presented here were collected from various sources – textbooks, articles, manuscripts, internet and newsletters. All the researchers and authors of the above mentioned sources are greatly acknowledged.