

CLA 3 - PRACTICAL COMPONENT - MCQ

pp0783@srmist.edu.in [Switch account](#)



Your email will be recorded when you submit this form

* Required

MCQ - Answer all the questions

The electrolyte solution within the glass electrode is *

- ☐ saturated KCl
- ☐ dilute KCl
- ☐ dilute NaCl
- ☒ saturated NaCl

A buffer solution comprises which one of the following? *

- ☐ a weak acid in solution
- ☒ a weak acid and its conjugate base in solution
- ☐ a strong acid and its conjugate base in solution
- ☐ a weak base in solution



In conductometric titration the titrant concentration must be as that of the titrand *

- ☒ same
- ☐ 5 times
- ☐ 10 times
- ☐ does not depend on the concentration

Soft water + Buffer + EBT gives *

- ☒ Appearance of wine-red colour
- ☐ Formation of weak complex
- ☐ Formation of reddish brown precipitate
- ☐ Appearance of steel blue colour

Titre value of methyl orange corresponds to *

- ☐ sodium hydroxide plus half the carbonate
- ☐ sodium hydroxide only
- ☐ sodium hydroxide and sodium carbonate
- ☒ sodium carbonate only



On increasing the temperature, the viscosity of the fluid _____ *

- ☐ Initially decreases then increases
- ☒ Increases
- ☐ Decreases
- ☐ No change is observed

When pH is below 8.5 the indicator ---- is colourless *

- ☐ Methyl orange
- ☐ Eriochrome black - T
- ☒ Phenolphthalein
- ☐ Potassium chromate

Which of the following chemical agent is added during the estimation of Fe(II) ions by potentiometry to avoid the hydrolysis reaction during the titration? *

- ☐ dilute HCl
- ☐ Conc. HCl
- ☒ dilute sulphuric acid
- ☐ dilute nitric acid



End point sharpness in precipitation titration will be improved if *

- ☐ reagent concentration decreases and product precipitate K_{sp} decreases
- ☒ reagent concentration increases and product precipitate K_{sp} increases
- ☐ reagent concentration decreases and product precipitate K_{sp} increases
- ☐ reagent concentration increases and product precipitate K_{sp} decreases

At the same concentration and temperature, dilute aqueous solution of strong acid will conduct electricity.... *

- ☒ better than dilute aqueous solution of weak acid
- ☐ as much as dilute aqueous solution of weak acid
- ☐ lower than the dilute aqueous solution of weak acid
- ☐ two-fold higher than the weak acid

All of the following statements are correct regarding potentiometric titration except *

- ☐ They are suitable for colored or turbid solutions
- ☐ The results obtained are accurate
- ☐ Acid base titration can also be carried out by potentiometry
- ☒ The EMF of the cell is zero at the equivalence point



Which of the following is not a primary standard? *

- ☐ Silver nitrate
- ☒ Oxalic acid
- ☐ Sodium chloride
- ☐ Anhydrous sodium carbonate

Glass electrode do not permit the direct potentiometric measurement of ion *

- ☐ Na⁺
- ☐ H⁺
- ☐ Ca²⁺
- ☒ None of the above

In determination of mixture of bases by titration method, the amount of Sodium Hydroxide is calculated as---. *

- ☐ $N \times \text{Equivalent mass of Sodium Carbonate} / 10$
- ☐ $N [\text{OH and CO}_3^{2-} \text{ portion}] \times \text{Equivalent mass of Sodium Hydroxide and Sodium carbonate} / 10$
- ☐ $N [\text{OH portion}] \times \text{Equivalent mass of Sodium Hydroxide} / 10$
- ☒ $N [\text{CO}_3^{2-} \text{ portion}] \times \text{Equivalent mass of Sodium carbonate} / 10$



25ml of 0.1 N NaOH solution is required for complete neutralisation of 0.32 g of an acid, What will be the molecular weight of the acid ? *

- ☐ 64
- ☒ 128
- ☐ 32
- ☐ 16

To prepare 25 ml of 0.2 % diluted solution from a 1% solution, we need *

- ☐ 10 ml of 1 % solution
- ☒ 5 ml of 1 % solution
- ☐ 2.5 ml of 1 % solution
- ☐ 7.5 ml of 1 % solution

Which of the following represents the equivalence point in the graph of EMF vs volume of titrant? *

- ☒ Point at the highest EMF
- ☐ Point at the lowest EMF
- ☐ Point at the greatest magnitude of the slope of the curve
- ☐ Point at the least magnitude of the slope of the curve



Name the reference electrode and working electrode used in the estimation of Fe(II) ions by potentiometry *

- ☐ Platinum electrode and Standard Calomel Electrode
- ☒ Standard Calomel Electrode and Platinum electrode
- ☐ Standard Calomel Electrode and Glass electrode
- ☐ Glass electrode and Platinum electrode

If 20 g of NaOH is dissolved in 500 mL of distilled water, then what is the concentration of the solution? *

- ☒ 1 N
- ☐ 0.05 N
- ☐ 0.5 N
- ☐ 0.025 N

The titration of FAS vs Potassium dichromate yields the following : $M\text{Fe}^{2+} + V\text{Fe}^{2+} = X M\text{Cr}_2\text{O}_7 + V\text{Cr}_2\text{O}_7$. What is X ? *

- ☐ 3
- ☒ 1
- ☐ 6
- ☐ 4



Why do we have to standardize AgNO_3 solution? *

- ☐ To calculate the normality of AgCl
- ☐ To find the normality of NaCl
- ☐ To calculate the volume of NaCl
- ☒ To find the normality of AgNO_3

In EDTA method, the formation constant of the metal-EDTA should be *

- ☒ above 10^4
- ☐ below 10^4
- ☐ at equilibrium = 0
- ☐ not dependent on formation constant

What is the pH range in which chloride determination using Mohr's method is conducted? *

- ☐ less than 4
- ☐ neutral pH
- ☐ 9-12
- ☒ 6-9



The significance of first derivative and second derivative plot in potentiometric titration is *

- ☐ To get additional information about the redox reaction
- ☒ To get more accurate equivalence point in case of colored and dilute solutions
- ☐ To get the value of standard electrode potential
- ☐ To get the voltage of reference electrode

In the plot of reduced viscosity versus concentration, the Y-intercept is *

- ☐ specific viscosity
- ☒ intrinsic viscosity
- ☐ inherent viscosity
- ☐ relative viscosity

Page 3 of 3

[Back](#)

[Submit](#)

[Clear form](#)

Never submit passwords through Google Forms.

This form was created inside of SRM Institute of Science and Technology. [Report Abuse](#)

Google Forms

